

2018 PRODUCT LIST

Columns and Accessories



NEW CHIRALPAK® IG-U

CHIRALPAK IH AVAILABLE
JANUARY 2018

CHIRALPAK® IG-U

Amylose tris
(3-chloro-5-methylphenylcarbamate)
Part No. 87U83
Particle size: 1.6 μm
Internal diameter: 3.0 mm
Length: 100 mm



CHIRAL
TECHNOLOGIES EUROPE
DAICEL GROUP

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (µ)	TYPE
-----------	-------	----------------	------------	----------------------	------

IMMOBILISED POLYSACCHARIDE CSPs AND COLUMNS*

CHIRALPAK® IA

80020	CHIRALPAK® IA	-	Kg	20	CSP
80021	CHIRALPAK® IA	-	100 g	20	CSP
80224	CHIRALPAK® IA	150	4.6	20	Analytical
80225	CHIRALPAK® IA	250	4.6	20	Analytical
80256	CHIRALPAK® IA	500	50	20	Preparative
80266	CHIRALPAK® IA	500	100	20	Preparative
8000E	CHIRALPAK® IA-10	-	Kg	10	CSP
80G24	CHIRALPAK® IA-10 SFC	150	4.6	10	SFC Analytical
80G25	CHIRALPAK® IA-10 SFC	250	4.6	10	SFC Analytical
80G35	CHIRALPAK® IA-10 SFC	250	10	10	SFC Semi-Prep
80G45	CHIRALPAK® IA-10 SFC	250	20	10	SFC Semi-Prep
80G75	CHIRALPAK® IA-10 SFC	250	30	10	SFC Semi-Prep
80311	CHIRALPAK® IA	10	4	5	Guard Cartridges (x 3)
80321	CHIRALPAK® IA	50	4.6	5	Analytical
80323	CHIRALPAK® IA	100	4.6	5	Analytical
80324	CHIRALPAK® IA	150	4.6	5	Analytical
80325	CHIRALPAK® IA	250	4.6	5	Analytical
80335	CHIRALPAK® IA	250	10	5	Semi-Prep
80337	CHIRALPAK® IA	20	10	5	Guard Semi-Prep
80345	CHIRALPAK® IA	250	20	5	Semi-Prep
80355	CHIRALPAK® IA	250	50	5	Preparative
80375	CHIRALPAK® IA	250	30	5	Semi-Prep
80394	CHIRALPAK® IA	150	2.1	5	Analytical
80423	CHIRALPAK® IA SFC	100	4.6	5	SFC Analytical
80435	CHIRALPAK® IA SFC	250	10	5	SFC Semi-Prep
80445	CHIRALPAK® IA SFC	250	20	5	SFC Semi-Prep
80455	CHIRALPAK® IA SFC	250	50	5	SFC Preparative
80475	CHIRALPAK® IA SFC	250	30	5	SFC Semi-Prep

CHIRALPAK® IA-3

80511	CHIRALPAK® IA-3	10	4	3	Guard Cartridges (x 3)
80522	CHIRALPAK® IA-3	50	4.6	3	Analytical
80523	CHIRALPAK® IA-3	100	4.6	3	Analytical
80524	CHIRALPAK® IA-3	150	4.6	3	Analytical
80525	CHIRALPAK® IA-3	250	4.6	3	Analytical
80592	CHIRALPAK® IA-3	50	2.1	3	Analytical
80593	CHIRALPAK® IA-3	100	2.1	3	Analytical
80594	CHIRALPAK® IA-3	150	2.1	3	Analytical
80595	CHIRALPAK® IA-3	250	2.1	3	Analytical
80S82	CHIRALPAK® IA-3 SFC	50	3	3	SFC Analytical
80S83	CHIRALPAK® IA-3 SFC	100	3	3	SFC Analytical
80S84	CHIRALPAK® IA-3 SFC	150	3	3	SFC Analytical

CHIRALPAK® IA-U

80U82	CHIRALPAK® IA-U	50	3	1.6	Analytical
80U83	CHIRALPAK® IA-U	100	3	1.6	Analytical

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (µ)	TYPE
-----------	-------	----------------	------------	----------------------	------

CHIRALPAK® IB

81311	CHIRALPAK® IB	10	4	5	Guard Cartridges (x 3)
81322	CHIRALPAK® IB	50	4.6	5	Analytical
81323	CHIRALPAK® IB	100	4.6	5	Analytical
81324	CHIRALPAK® IB	150	4.6	5	Analytical
81325	CHIRALPAK® IB	250	4.6	5	Analytical
81335	CHIRALPAK® IB	250	10	5	Semi-Prep
81337	CHIRALPAK® IB	20	10	5	Guard Semi-Prep
81345	CHIRALPAK® IB	250	20	5	Semi-Prep
81355	CHIRALPAK® IB	250	50	5	Preparative
81375	CHIRALPAK® IB	250	30	5	Semi-Prep
81394	CHIRALPAK® IB	150	2.1	5	Analytical
81423	CHIRALPAK® IB SFC	100	4.6	5	SFC Analytical
81435	CHIRALPAK® IB SFC	250	10	5	SFC Semi-Prep
81445	CHIRALPAK® IB SFC	250	20	5	SFC Semi-Prep
81455	CHIRALPAK® IB SFC	250	50	5	SFC Preparative
81475	CHIRALPAK® IB SFC	250	30	5	SFC Semi-Prep

CHIRALPAK® IB-3

81511	CHIRALPAK® IB-3	10	4	3	Guard Cartridges (x 3)
81522	CHIRALPAK® IB-3	50	4.6	3	Analytical
81523	CHIRALPAK® IB-3	100	4.6	3	Analytical
81524	CHIRALPAK® IB-3	150	4.6	3	Analytical
81525	CHIRALPAK® IB-3	250	4.6	3	Analytical
81592	CHIRALPAK® IB-3	50	2.1	3	Analytical
81593	CHIRALPAK® IB-3	100	2.1	3	Analytical
81594	CHIRALPAK® IB-3	150	2.1	3	Analytical
81595	CHIRALPAK® IB-3	250	2.1	3	Analytical
81S82	CHIRALPAK® IB-3 SFC	50	3	3	SFC Analytical
81S83	CHIRALPAK® IB-3 SFC	100	3	3	SFC Analytical
81S84	CHIRALPAK® IB-3 SFC	150	3	3	SFC Analytical

CHIRALPAK® IB-U

81U82	CHIRALPAK® IB-U	50	3.0	1.6	Analytical
81U83	CHIRALPAK® IB-U	100	3.0	1.6	Analytical

CHIRALPAK® IB N-20

88020	CHIRALPAK® IB N-20	-	Kg	20	CSP
88021	CHIRALPAK® IB N-20	-	100 g	20	CSP
88224	CHIRALPAK® IB N-20	150	4.6	20	Analytical
88225	CHIRALPAK® IB N-20	250	4.6	20	Analytical
88256	CHIRALPAK® IB N-20	500	50	20	Preparative
88266	CHIRALPAK® IB N-20	500	100	20	Preparative

CHIRALPAK® IB N-10

8800E	CHIRALPAK® IB N-10	-	Kg	10	CSP
88G24	CHIRALPAK® IB N-10 SFC	150	4.6	10	SFC Analytical
88G25	CHIRALPAK® IB N-10 SFC	250	4.6	10	SFC Analytical
88G35	CHIRALPAK® IB N-10 SFC	250	10	10	SFC Semi-Prep
88G45	CHIRALPAK® IB N-10 SFC	250	21	10	SFC Semi-Prep
88G75	CHIRALPAK® IB N-10 SFC	250	30	10	SFC Semi-Prep

*Immobilised polysaccharide columns are compatible for use in either normal or reverse phase chromatographic conditions.

2018 Product List

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
CHIRALPAK® IB N-5					
88311	CHIRALPAK® IB N-5	10	4	5	Guard Cartridges (x 3)
88324	CHIRALPAK® IB N-5	150	4.6	5	Analytical
88325	CHIRALPAK® IB N-5	250	4.6	5	Analytical
88335	CHIRALPAK® IB N-5	250	10	5	Semi-Prep
88337	CHIRALPAK® IB N-5	20	10	5	Guard Semi-Prep
88345	CHIRALPAK® IB N-5	250	20	5	Semi-Prep
88375	CHIRALPAK® IB N-5	250	30	5	Semi-Prep
88394	CHIRALPAK® IB N-5	150	2.1	5	Analytical
88422	CHIRALPAK® IB N-5 SFC	50	4.6	5	SFC Analytical
88423	CHIRALPAK® IB N-5 SFC	100	4.6	5	SFC Analytical
88435	CHIRALPAK® IB N-5 SFC	250	10	5	SFC Semi-Prep
88445	CHIRALPAK® IB N-5 SFC	250	21	5	SFC Semi-Prep
88455	CHIRALPAK® IB N-5 SFC	250	50	5	SFC Preparative
88475	CHIRALPAK® IB N-5 SFC	250	30	5	SFC Semi-Prep

CHIRALPAK® IB N-3

88511	CHIRALPAK® IB N-3	10	4	3	Guard Cartridges (x 3)
88522	CHIRALPAK® IB N-3	50	4.6	3	Analytical
88523	CHIRALPAK® IB N-3	100	4.6	3	Analytical
88524	CHIRALPAK® IB N-3	150	4.6	3	Analytical
88525	CHIRALPAK® IB N-3	250	4.6	3	Analytical
88592	CHIRALPAK® IB N-3	50	2.1	3	Analytical
88593	CHIRALPAK® IB N-3	100	2.1	3	Analytical
88594	CHIRALPAK® IB N-3	150	2.1	3	Analytical
88595	CHIRALPAK® IB N-3	250	2.1	3	Analytical
88S82	CHIRALPAK® IB N-3 SFC	50	3	3	SFC Analytical
88S83	CHIRALPAK® IB N-3 SFC	100	3	3	SFC Analytical
88S84	CHIRALPAK® IB N-3 SFC	150	3	3	SFC Analytical

CHIRALCEL® OD-I

82020	CHIRALCEL® OD-I	-	Kg	20	CSP
82021	CHIRALCEL® OD-I	-	100 g	20	CSP
82224	CHIRALCEL® OD-I	150	4.6	20	Analytical
82225	CHIRALCEL® OD-I	250	4.6	20	Analytical
82256	CHIRALCEL® OD-I	500	50	20	Preparative
82266	CHIRALCEL® OD-I	500	100	20	Preparative

CHIRALPAK® IC

83020	CHIRALPAK® IC	-	Kg	20	CSP
83021	CHIRALPAK® IC	-	100 g	20	CSP
83224	CHIRALPAK® IC	150	4.6	20	Analytical
83225	CHIRALPAK® IC	250	4.6	20	Analytical
83256	CHIRALPAK® IC	500	50	20	Preparative
83266	CHIRALPAK® IC	500	100	20	Preparative
8300E	CHIRALPAK® IC-10	-	Kg	10	CSP
83G24	CHIRALPAK® IC-10 SFC	150	4.6	10	SFC Analytical
83G25	CHIRALPAK® IC-10 SFC	250	4.6	10	SFC Analytical
83G35	CHIRALPAK® IC-10 SFC	250	10	10	SFC Semi-Prep
83G45	CHIRALPAK® IC-10 SFC	250	20	10	SFC Semi-Prep
83G75	CHIRALPAK® IC-10 SFC	250	30	10	SFC Semi-Prep
83311	CHIRALPAK® IC	10	4	5	Guard Cartridges (x 3)

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
83322	CHIRALPAK® IC	50	4.6	5	Analytical
83323	CHIRALPAK® IC	100	4.6	5	Analytical
83324	CHIRALPAK® IC	150	4.6	5	Analytical
83325	CHIRALPAK® IC	250	4.6	5	Analytical
83335	CHIRALPAK® IC	250	10	5	Semi-Prep
83337	CHIRALPAK® IC	20	10	5	Guard Semi-Prep
83345	CHIRALPAK® IC	250	20	5	Semi-Prep
83355	CHIRALPAK® IC	250	50	5	Preparative
83375	CHIRALPAK® IC	250	30	5	Semi-Prep
83394	CHIRALPAK® IC	150	2.1	5	Analytical
83423	CHIRALPAK® IC SFC	100	4.6	5	SFC Analytical
83435	CHIRALPAK® IC SFC	250	10	5	SFC Semi-Prep
83445	CHIRALPAK® IC SFC	250	20	5	SFC Semi-Prep
83455	CHIRALPAK® IC SFC	250	50	5	SFC Preparative
83475	CHIRALPAK® IC SFC	250	30	5	SFC Semi-Prep

CHIRALPAK® IC-3

83511	CHIRALPAK® IC-3	10	4	3	Guard Cartridges (x 3)
83522	CHIRALPAK® IC-3	50	4.6	3	Analytical
83523	CHIRALPAK® IC-3	100	4.6	3	Analytical
83524	CHIRALPAK® IC-3	150	4.6	3	Analytical
83525	CHIRALPAK® IC-3	250	4.6	3	Analytical
83592	CHIRALPAK® IC-3	50	2.1	3	Analytical
83593	CHIRALPAK® IC-3	100	2.1	3	Analytical
83594	CHIRALPAK® IC-3	150	2.1	3	Analytical
83595	CHIRALPAK® IC-3	250	2.1	3	Analytical
80S82	CHIRALPAK® IC-3 SFC	50	3	3	SFC Analytical
83S83	CHIRALPAK® IC-3 SFC	100	3	3	SFC Analytical
83S84	CHIRALPAK® IC-3 SFC	150	3	3	SFC Analytical

CHIRALPAK® IC-U

83U82	CHIRALPAK® IC-U	50	3	1.6	Analytical
83U83	CHIRALPAK® IC-U	100	3	1.6	Analytical

CHIRALPAK® ID

84020	CHIRALPAK® ID	-	Kg	20	CSP
84021	CHIRALPAK® ID	-	100 g	20	CSP
84224	CHIRALPAK® ID	150	4.6	20	Analytical
84225	CHIRALPAK® ID	250	4.6	20	Analytical
84256	CHIRALPAK® ID	500	50	20	Preparative
84266	CHIRALPAK® ID	500	100	20	Preparative
8400E	CHIRALPAK® ID-10	-	Kg	10	CSP
84G24	CHIRALPAK® ID-10 SFC	150	4.6	10	SFC Analytical
84G25	CHIRALPAK® ID-10 SFC	250	4.6	10	SFC Analytical
84G35	CHIRALPAK® ID-10 SFC	250	10	10	SFC Semi-Prep
84G45	CHIRALPAK® ID-10 SFC	250	20	10	SFC Semi-Prep
84G75	CHIRALPAK® ID-10 SFC	250	30	10	SFC Semi-Prep
84311	CHIRALPAK® ID	10	4	5	Guard Cartridges (x 3)
84321	CHIRALPAK® ID	50	4.6	5	Analytical
84323	CHIRALPAK® ID	100	4.6	5	Analytical
84324	CHIRALPAK® ID	150	4.6	5	Analytical
84325	CHIRALPAK® ID	250	4.6	5	Analytical

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
84335	CHIRALPAK® ID	250	10	5	Semi-Prep
84337	CHIRALPAK® ID	20	10	5	Guard Semi-Prep
84345	CHIRALPAK® ID	250	20	5	Semi-Prep
84355	CHIRALPAK® ID	250	50	5	Preparative
84375	CHIRALPAK® ID	250	30	5	Semi-Prep
84394	CHIRALPAK® ID	150	2.1	5	Analytical
84423	CHIRALPAK® ID SFC	100	4.6	5	SFC Analytical
84435	CHIRALPAK® ID SFC	250	10	5	SFC Semi-Prep
84445	CHIRALPAK® ID SFC	250	20	5	SFC Semi-Prep
84455	CHIRALPAK® ID SFC	250	50	5	SFC Preparative
84475	CHIRALPAK® ID SFC	250	30	5	SFC Semi-Prep

CHIRALPAK® ID-3

84511	CHIRALPAK® ID-3	10	4	3	Guard Cartridges (x 3)
84522	CHIRALPAK® ID-3	50	4.6	3	Analytical
84523	CHIRALPAK® ID-3	100	4.6	3	Analytical
84524	CHIRALPAK® ID-3	150	4.6	3	Analytical
84525	CHIRALPAK® ID-3	250	4.6	3	Analytical
84592	CHIRALPAK® ID-3	50	2.1	3	Analytical
84593	CHIRALPAK® ID-3	100	2.1	3	Analytical
84594	CHIRALPAK® ID-3	150	2.1	3	Analytical
84595	CHIRALPAK® ID-3	250	2.1	3	Analytical
84S82	CHIRALPAK® ID-3 SFC	50	3	3	SFC Analytical
84S83	CHIRALPAK® ID-3 SFC	100	3	3	SFC Analytical
84S84	CHIRALPAK® ID-3 SFC	150	3	3	SFC Analytical

CHIRALPAK® IE

85020	CHIRALPAK® IE	-	Kg	20	CSP
85021	CHIRALPAK® IE	-	100 g	20	CSP
85224	CHIRALPAK® IE	150	4.6	20	Analytical
85225	CHIRALPAK® IE	250	4.6	20	Analytical
85256	CHIRALPAK® IE	500	50	20	Preparative
85266	CHIRALPAK® IE	500	100	20	Preparative
8500E	CHIRALPAK® IE-10	-	Kg	10	CSP
85G24	CHIRALPAK® IE-10 SFC	150	4.6	10	SFC Analytical
85G25	CHIRALPAK® IE-10 SFC	250	4.6	10	SFC Analytical
85G35	CHIRALPAK® IE-10 SFC	250	10	10	SFC Semi-Prep
85G45	CHIRALPAK® IE-10 SFC	250	20	10	SFC Semi-Prep
85G75	CHIRALPAK® IE-10 SFC	250	30	10	SFC Semi-Prep
85311	CHIRALPAK® IE	10	4	5	Guard Cartridges (x 3)
85321	CHIRALPAK® IE	50	4.6	5	Analytical
85323	CHIRALPAK® IE	100	4.6	5	Analytical
85324	CHIRALPAK® IE	150	4.6	5	Analytical
85325	CHIRALPAK® IE	250	4.6	5	Analytical
85335	CHIRALPAK® IE	250	10	5	Semi-Prep
85337	CHIRALPAK® IE	20	10	5	Guard Semi-Prep
85345	CHIRALPAK® IE	250	20	5	Semi-Prep
85355	CHIRALPAK® IE	250	50	5	Preparative
85375	CHIRALPAK® IE	250	30	5	Semi-Prep
85394	CHIRALPAK® IE	150	2.1	5	Analytical
85423	CHIRALPAK® IE SFC	100	4.6	5	SFC Analytical
85435	CHIRALPAK® IE SFC	250	10	5	SFC Semi-Prep

85445	CHIRALPAK® IE SFC	250	20	5	SFC Semi-Prep
85455	CHIRALPAK® IE SFC	250	50	5	SFC Preparative
85475	CHIRALPAK® IE SFC	250	30	5	SFC Semi-Prep

CHIRALPAK® IE-3

85511	CHIRALPAK® IE-3	10	4	3	Guard Cartridges (x 3)
85522	CHIRALPAK® IE-3	50	4.6	3	Analytical
85523	CHIRALPAK® IE-3	100	4.6	3	Analytical
85524	CHIRALPAK® IE-3	150	4.6	3	Analytical
85525	CHIRALPAK® IE-3	250	4.6	3	Analytical
85592	CHIRALPAK® IE-3	50	2.1	3	Analytical
85593	CHIRALPAK® IE-3	100	2.1	3	Analytical
85594	CHIRALPAK® IE-3	150	2.1	3	Analytical
85595	CHIRALPAK® IE-3	250	2.1	3	Analytical
85S82	CHIRALPAK® IE-3 SFC	50	3	3	SFC Analytical
85S83	CHIRALPAK® IE-3 SFC	100	3	3	SFC Analytical
85S84	CHIRALPAK® IE-3 SFC	150	3	3	SFC Analytical

CHIRALPAK® IF

86020	CHIRALPAK® IF	-	Kg	20	CSP
86021	CHIRALPAK® IF	-	100 g	20	CSP
86224	CHIRALPAK® IF	150	4.6	20	Analytical
86225	CHIRALPAK® IF	250	4.6	20	Analytical
86256	CHIRALPAK® IF	500	50	20	Preparative
86266	CHIRALPAK® IF	500	100	20	Preparative
8600E	CHIRALPAK® IF-10	-	Kg	10	CSP
86G24	CHIRALPAK® IF-10 SFC	150	4.6	10	SFC Analytical
86G25	CHIRALPAK® IF-10 SFC	250	4.6	10	SFC Analytical
86G35	CHIRALPAK® IF-10 SFC	250	10	10	SFC Semi-Prep
86G45	CHIRALPAK® IF-10 SFC	250	20	10	SFC Semi-Prep
86G75	CHIRALPAK® IF-10 SFC	250	30	10	SFC Semi-Prep
86311	CHIRALPAK® IF	10	4	5	Guard Cartridges (x 3)
86321	CHIRALPAK® IF	50	4.6	5	Analytical
86323	CHIRALPAK® IF	100	4.6	5	Analytical
86324	CHIRALPAK® IF	150	4.6	5	Analytical
86325	CHIRALPAK® IF	250	4.6	5	Analytical
86335	CHIRALPAK® IF	250	10	5	Semi-Prep
86337	CHIRALPAK® IF	20	10	5	Guard Semi-Prep
86345	CHIRALPAK® IF	250	20	5	Semi-Prep
86355	CHIRALPAK® IF	250	50	5	Preparative
86375	CHIRALPAK® IF	250	30	5	Semi-Prep
86394	CHIRALPAK® IF	150	2.1	5	Analytical
86423	CHIRALPAK® IF SFC	100	4.6	5	SFC Analytical
86435	CHIRALPAK® IF SFC	250	10	5	SFC Semi-Prep
86445	CHIRALPAK® IF SFC	250	20	5	SFC Semi-Prep
86455	CHIRALPAK® IF SFC	250	50	5	SFC Preparative
86475	CHIRALPAK® IF SFC	250	30	5	SFC Semi-Prep

CHIRALPAK® IF-3

86511	CHIRALPAK® IF-3	10	4	3	Guard Cartridges (x 3)
86522	CHIRALPAK® IF-3	50	4.6	3	Analytical
86523	CHIRALPAK® IF-3	100	4.6	3	Analytical

2018 Product List

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
86524	CHIRALPAK® IF-3	150	4.6	3	Analytical
86525	CHIRALPAK® IF-3	250	4.6	3	Analytical
86592	CHIRALPAK® IF-3	50	2.1	3	Analytical
86593	CHIRALPAK® IF-3	100	2.1	3	Analytical
86594	CHIRALPAK® IF-3	150	2.1	3	Analytical
86595	CHIRALPAK® IF-3	250	2.1	3	Analytical
86S82	CHIRALPAK® IF-3 SFC	50	3	3	SFC Analytical
86S83	CHIRALPAK® IF-3 SFC	100	3	3	SFC Analytical
86S84	CHIRALPAK® IF-3 SFC	150	3	3	SFC Analytical

CHIRALPAK® IG

87020	CHIRALPAK® IG	-	Kg	20	CSP
87021	CHIRALPAK® IG	-	100 g	20	CSP
87224	CHIRALPAK® IG	150	4.6	20	Analytical
87225	CHIRALPAK® IG	250	4.6	20	Analytical
87256	CHIRALPAK® IG	500	50	20	Preparative
87266	CHIRALPAK® IG	500	100	20	Preparative
8700E	CHIRALPAK® IG-10	-	Kg	10	CSP
87G24	CHIRALPAK® IG-10 SFC	150	4.6	10	SFC Analytical
87G25	CHIRALPAK® IG-10 SFC	250	4.6	10	SFC Analytical
87G35	CHIRALPAK® IG-10 SFC	250	10	10	SFC Semi-Prep
87G45	CHIRALPAK® IG-10 SFC	250	20	10	SFC Semi-Prep
87G75	CHIRALPAK® IG-10 SFC	250	30	10	SFC Semi-Prep
87311	CHIRALPAK® IG	10	4	5	Guard Cartridges (x 3)
87321	CHIRALPAK® IG	50	4.6	5	Analytical
87323	CHIRALPAK® IG	100	4.6	5	Analytical
87324	CHIRALPAK® IG	150	4.6	5	Analytical
87325	CHIRALPAK® IG	250	4.6	5	Analytical
87335	CHIRALPAK® IG	250	10	5	Semi-Prep
87337	CHIRALPAK® IG	20	10	5	Guard Semi-Prep
87345	CHIRALPAK® IG	250	20	5	Semi-Prep
87355	CHIRALPAK® IG	250	50	5	Preparative
87375	CHIRALPAK® IG	250	30	5	Semi-Prep
87394	CHIRALPAK® IG	150	2.1	5	Analytical
87423	CHIRALPAK® IG SFC	100	4.6	5	SFC Analytical
87435	CHIRALPAK® IG SFC	250	10	5	SFC Semi-Prep
87445	CHIRALPAK® IG SFC	250	20	5	SFC Semi-Prep
87455	CHIRALPAK® IG SFC	250	50	5	SFC Preparative
87475	CHIRALPAK® IG SFC	250	30	5	SFC Semi-Prep

CHIRALPAK® IG-3

87511	CHIRALPAK® IG-3	10	4	3	Guard Cartridges (x 3)
87522	CHIRALPAK® IG-3	50	4.6	3	Analytical
87523	CHIRALPAK® IG-3	100	4.6	3	Analytical
87524	CHIRALPAK® IG-3	150	4.6	3	Analytical
87525	CHIRALPAK® IG-3	250	4.6	3	Analytical
87592	CHIRALPAK® IG-3	50	2.1	3	Analytical
87593	CHIRALPAK® IG-3	100	2.1	3	Analytical
87594	CHIRALPAK® IG-3	150	2.1	3	Analytical
87595	CHIRALPAK® IG-3	250	2.1	3	Analytical
87S82	CHIRALPAK® IG-3 SFC	50	3	3	SFC Analytical
87S83	CHIRALPAK® IG-3 SFC	100	3	3	SFC Analytical
87S84	CHIRALPAK® IG-3 SFC	150	3	3	SFC Analytical

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
CHIRALPAK® IG-U					
87U82	CHIRALPAK® IG-U	50	3	1.6	Analytical
87U83	CHIRALPAK® IG-U	100	3	1.6	Analytical

CHIRALPAK® IH

89311	CHIRALPAK® IH	10	4	5	Guard Cartridges (x 3)
89321	CHIRALPAK® IH	50	4.6	5	Analytical
89323	CHIRALPAK® IH	100	4.6	5	Analytical
89324	CHIRALPAK® IH	150	4.6	5	Analytical
89325	CHIRALPAK® IH	250	4.6	5	Analytical
89335	CHIRALPAK® IH	250	10	5	Semi-Prep
89337	CHIRALPAK® IH	20	10	5	Guard Semi-Prep
89345	CHIRALPAK® IH	250	20	5	Semi-Prep
89355	CHIRALPAK® IH	250	50	5	Preparative
89375	CHIRALPAK® IH	250	30	5	Semi-Prep
89394	CHIRALPAK® IH	150	2.1	5	Analytical
89423	CHIRALPAK® IH SFC	100	4.6	5	SFC Analytical
89435	CHIRALPAK® IH SFC	250	10	5	SFC Semi-Prep
89445	CHIRALPAK® IH SFC	250	20	5	SFC Semi-Prep
89455	CHIRALPAK® IH SFC	250	50	5	SFC Preparative
89475	CHIRALPAK® IH SFC	250	30	5	SFC Semi-Prep

CHIRALCEL® OA

11022	CHIRALCEL® OA	50	4.6	10	Guard
11025	CHIRALCEL® OA	250	4.6	10	HPLC Analytical

CHIRALCEL® OB

12022	CHIRALCEL® OB	50	4.6	10	Guard
12025	CHIRALCEL® OB	250	4.6	10	HPLC Analytical

CHIRALCEL® OB-H

12311	CHIRALCEL® OB-H	10	4	5	Guard Cartridges (x 3)
12324	CHIRALCEL® OB-H	150	4.6	5	Analytical
12325	CHIRALCEL® OB-H	250	4.6	5	Analytical

CHIRALCEL® OC

13022	CHIRALCEL® OC	50	4.6	10	Guard
13025	CHIRALCEL® OC	250	4.6	10	HPLC Analytical

CHIRALCEL® OC-H

13325	CHIRALCEL® OC-H	250	4.6	5	Analytical
-------	-----------------	-----	-----	---	------------

CHIRALCEL® OD

14020	CHIRALCEL® OD	-	Kg	20	CSP
14021	CHIRALCEL® OD	-	100 g	20	CSP
14024	CHIRALCEL® OD	150	4.6	10	HPLC Analytical
14025	CHIRALCEL® OD	250	4.6	10	HPLC Analytical
14224	CHIRALCEL® OD	150	4.6	20	Analytical
14225	CHIRALCEL® OD	250	4.6	20	Analytical
14256	CHIRALCEL® OD	500	50	20	Preparative
14266	CHIRALCEL® OD	500	100	20	Preparative
14032	CHIRALCEL® OD	50	10	10	Guard Semi-Prep
14035	CHIRALCEL® OD	250	10	10	Semi-Prep
14042	CHIRALCEL® OD	50	20	10	Guard Semi-Prep

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
1400E	CHIRALCEL® OD-10	-	Kg	10	CSP
14G22	CHIRALCEL® OD-10 SFC	50	4.6	10	SFC Analytical
14G24	CHIRALCEL® OD-10 SFC	150	4.6	10	SFC Analytical
14G25	CHIRALCEL® OD-10 SFC	250	4.6	10	SFC Analytical
14G94	CHIRALCEL® OD-10 SFC	150	2.1	10	SFC Analytical
14G35	CHIRALCEL® OD-10 SFC	250	10	10	SFC Semi-Prep
14G45	CHIRALCEL® OD-10 SFC	250	20	10	SFC Semi-Prep
14G55	CHIRALCEL® OD-10 SFC	250	50	10	SFC Preparative
14G75	CHIRALCEL® OD-10 SFC	250	30	10	SFC Semi-Prep

CHIRALCEL® OD-H

14311	CHIRALCEL® OD-H	10	4	5	Guard Cartridges (x 3)
14321	CHIRALCEL® OD-H	50	4.6	5	Analytical
14323	CHIRALCEL® OD-H	100	4.6	5	Analytical
14324	CHIRALCEL® OD-H	150	4.6	5	Analytical
14325	CHIRALCEL® OD-H	250	4.6	5	Analytical
14335	CHIRALCEL® OD-H	250	10	5	Semi-Prep
14337	CHIRALCEL® OD-H	20	10	5	Guard Semi-Prep
14345	CHIRALCEL® OD-H	250	20	5	Semi-Prep
14355	CHIRALCEL® OD-H	250	50	5	Preparative
14375	CHIRALCEL® OD-H	250	30	5	Semi-Prep
14394	CHIRALCEL® OD-H	150	2.1	5	Analytical
14423	CHIRALCEL® OD-H SFC	100	4.6	5	SFC Analytical
14435	CHIRALCEL® OD-H SFC	250	10	5	SFC Semi-Prep
14445	CHIRALCEL® OD-H SFC	250	20	5	SFC Semi-Prep
14455	CHIRALCEL® OD-H SFC	250	50	5	SFC Preparative
14475	CHIRALCEL® OD-H SFC	250	30	5	SFC Semi-Prep

CHIRALCEL® OD-3

14511	CHIRALCEL® OD-3	10	4	3	Guard Cartridges (x 3)
14522	CHIRALCEL® OD-3	50	4.6	3	Analytical
14523	CHIRALCEL® OD-3	100	4.6	3	Analytical
14524	CHIRALCEL® OD-3	150	4.6	3	Analytical
14525	CHIRALCEL® OD-3	250	4.6	3	Analytical
14592	CHIRALCEL® OD-3	50	2.1	3	Analytical
14593	CHIRALCEL® OD-3	100	2.1	3	Analytical
14594	CHIRALCEL® OD-3	150	2.1	3	Analytical
14595	CHIRALCEL® OD-3	250	2.1	3	Analytical
14S82	CHIRALCEL® OD-3 SFC	50	3	3	SFC Analytical
14S83	CHIRALCEL® OD-3 SFC	100	3	3	SFC Analytical
14S84	CHIRALCEL® OD-3 SFC	150	3	3	SFC Analytical

CHIRALCEL® OD-R

14611	CHIRALCEL® OD-R	10	4	10	Guard Cartridges (x 3)
14625	CHIRALCEL® OD-R	250	4.6	10	HPLC Analytical

CHIRALCEL® OD-RH

14711	CHIRALCEL® OD-RH	10	4	5	Guard Cartridges (x 3)
14724	CHIRALCEL® OD-RH	150	4.6	5	Analytical
14794	CHIRALCEL® OD-RH	150	2.1	5	Analytical

CHIRALCEL® OD-3R

14811	CHIRALCEL® OD-3R	10	4	3	Guard Cartridges (x 3)
14822	CHIRALCEL® OD-3R	50	4.6	3	Analytical

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
14823	CHIRALCEL® OD-3R	100	4.6	3	Analytical
14824	CHIRALCEL® OD-3R	150	4.6	3	Analytical
14892	CHIRALCEL® OD-3R	50	2.1	3	Analytical
14893	CHIRALCEL® OD-3R	100	2.1	3	Analytical
14894	CHIRALCEL® OD-3R	150	2.1	3	Analytical
14895	CHIRALCEL® OD-3R	250	2.1	3	Analytical

CHIRALCEL® OF

15022	CHIRALCEL® OF	50	4.6	10	Guard
15025	CHIRALCEL® OF	250	4.6	10	HPLC Analytical

CHIRALCEL® OG

16022	CHIRALCEL® OG	50	4.6	10	Guard
16025	CHIRALCEL® OG	250	4.6	10	HPLC Analytical

CHIRALCEL® OJ

17020	CHIRALCEL® OJ	-	Kg	20	CSP
17021	CHIRALCEL® OJ	-	100 g	20	CSP
17224	CHIRALCEL® OJ	150	4.6	20	Analytical
17225	CHIRALCEL® OJ	250	4.6	20	Analytical
17256	CHIRALCEL® OJ	500	50	20	Preparative
17266	CHIRALCEL® OJ	500	100	20	Preparative
17024	CHIRALCEL® OJ	150	4.6	10	HPLC Analytical
17025	CHIRALCEL® OJ	250	4.6	10	HPLC Analytical
17032	CHIRALCEL® OJ	50	10	10	Guard Semi-Prep
17035	CHIRALCEL® OJ	250	10	10	Semi-Prep
17042	CHIRALCEL® OJ	50	20	10	Guard Semi-Prep
17045	CHIRALCEL® OJ	250	20	10	Semi-Prep
1700E	CHIRALCEL® OJ-10	-	Kg	10	CSP
17G22	CHIRALCEL® OJ-10 SFC	50	4.6	10	SFC Analytical
17G24	CHIRALCEL® OJ-10 SFC	150	4.6	10	SFC Analytical
17G25	CHIRALCEL® OJ-10 SFC	250	4.6	10	SFC Analytical
17G94	CHIRALCEL® OJ-10 SFC	150	2.1	10	SFC Analytical
17G35	CHIRALCEL® OJ-10 SFC	250	10	10	SFC Semi-Prep
17G45	CHIRALCEL® OJ-10 SFC	250	21	10	SFC Semi-Prep
17G55	CHIRALCEL® OJ-10 SFC	250	50	10	SFC Preparative
17G75	CHIRALCEL® OJ-10 SFC	250	30	10	SFC Semi-Prep

CHIRALCEL® OJ-H

17311	CHIRALCEL® OJ-H	10	4	5	Guard Cartridges (x 3)
17321	CHIRALCEL® OJ-H	50	4.6	5	Analytical
17323	CHIRALCEL® OJ-H	100	4.6	5	Analytical
17324	CHIRALCEL® OJ-H	150	4.6	5	Analytical
17325	CHIRALCEL® OJ-H	250	4.6	5	Analytical
17335	CHIRALCEL® OJ-H	250	10	5	Semi-Prep
17337	CHIRALCEL® OJ-H	20	10	5	Guard Semi-Prep
17345	CHIRALCEL® OJ-H	250	20	5	Semi-Prep
17355	CHIRALCEL® OJ-H	250	50	5	Preparative
17375	CHIRALCEL® OJ-H	250	30	5	Semi-Prep
17394	CHIRALCEL® OJ-H	150	2.1	5	Analytical
17423	CHIRALCEL® OJ-H SFC	100	4.6	5	SFC Analytical
17435	CHIRALCEL® OJ-H SFC	250	10	5	SFC Semi-Prep
17445	CHIRALCEL® OJ-H SFC	250	20	5	SFC Semi-Prep

2018 Product List

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
17455	CHIRALCEL® OJ-H SFC	250	50	5	SFC Preparative
17475	CHIRALCEL® OJ-H SFC	250	30	5	SFC Semi-Prep

CHIRALCEL® OJ-3

17511	CHIRALCEL® OJ-3	10	4	3	Guard Cartridges (x 3)
17522	CHIRALCEL® OJ-3	50	4.6	3	Analytical
17523	CHIRALCEL® OJ-3	100	4.6	3	Analytical
17524	CHIRALCEL® OJ-3	150	4.6	3	Analytical
17525	CHIRALCEL® OJ-3	250	4.6	3	Analytical
17592	CHIRALCEL® OJ-3	50	2.1	3	Analytical
17593	CHIRALCEL® OJ-3	100	2.1	3	Analytical
17594	CHIRALCEL® OJ-3	150	2.1	3	Analytical
17595	CHIRALCEL® OJ-3	250	2.1	3	Analytical
17S82	CHIRALCEL® OJ-3 SFC	50	3	3	SFC Analytical
17S83	CHIRALCEL® OJ-3 SFC	100	3	3	SFC Analytical
17S84	CHIRALCEL® OJ-3 SFC	150	3	3	SFC Analytical

CHIRALCEL® OJ-RH

17711	CHIRALCEL® OJ-RH	10	4	5	Guard Cartridges (x 3)
17724	CHIRALCEL® OJ-RH	150	4.6	5	Analytical
17794	CHIRALCEL® OJ-RH	150	2.1	5	Analytical

CHIRALCEL® OJ-3R

17811	CHIRALCEL® OJ-3R	10	4	3	Guard Cartridges (x 3)
17822	CHIRALCEL® OJ-3R	50	4.6	3	Analytical
17823	CHIRALCEL® OJ-3R	100	4.6	3	Analytical
17824	CHIRALCEL® OJ-3R	150	4.6	3	Analytical
17892	CHIRALCEL® OJ-3R	50	2.1	3	Analytical
17893	CHIRALCEL® OJ-3R	100	2.1	3	Analytical
17894	CHIRALCEL® OJ-3R	150	2.1	3	Analytical
17895	CHIRALCEL® OJ-3R	250	2.1	3	Analytical

CHIRALCEL® OK

18022	CHIRALCEL® OK	50	4.6	10	Guard
18025	CHIRALCEL® OK	250	4.6	10	HPLC Analytical

CHIRALCEL® OX

63020	CHIRALCEL® OX	-	1 Kg	20	CSP
63021	CHIRALCEL® OX	-	100	20	CSP
63224	CHIRALCEL® OX	150	4.6	20	Analytical
63225	CHIRALCEL® OX	250	4.6	20	Analytical
63226	CHIRALCEL® OX	300	4.6	20	Analytical
63235	CHIRALCEL® OX	250	10	20	Semi-Prep
63242	CHIRALCEL® OX	50	21	20	Semi-Prep Guard
63245	CHIRALCEL® OX	250	21	20	Semi-Prep
63255	CHIRALCEL® OX	250	50	20	Preparative
63256	CHIRALCEL® OX	500	50	20	Preparative
63266	CHIRALCEL® OX	500	100	20	Preparative
63275	CHIRALCEL® OX	250	30	20	Semi-Prep
63935	CHIRALCEL® OX SFC	250	10	20	SFC Semi-Prep
63945	CHIRALCEL® OX SFC	250	21	20	SFC Semi-Prep
63955	CHIRALCEL® OX SFC	250	50	20	SFC Preparative

63975	CHIRALCEL® OX SFC	250	30	20	SFC Semi-Prep
63G22	CHIRALCEL® OX-10 SFC	50	4.6	10	SFC Analytical
63G24	CHIRALCEL® OX-10 SFC	150	4.6	10	SFC Analytical
63G25	CHIRALCEL® OX-10 SFC	250	4.6	10	SFC Analytical
63G35	CHIRALCEL® OX-10 SFC	250	10	10	SFC Semi-Prep
63G45	CHIRALCEL® OX-10 SFC	250	20	10	SFC Semi-Prep
63G75	CHIRALCEL® OX-10 SFC	250	30	10	SFC Semi-Prep
63G55	CHIRALCEL® OX-10 SFC	250	50	10	SFC Preparative
63G94	CHIRALCEL® OX-10 SFC	150	2.1	10	SFC Analytical

CHIRALCEL® OX-H

63311	CHIRALCEL® OX-H	10	4	5	Guard Cartridges (x 3)
63321	CHIRALCEL® OX-H	50	4.6	5	Analytical
63323	CHIRALCEL® OX-H	100	4.6	5	Analytical
63324	CHIRALCEL® OX-H	150	4.6	5	Analytical
63325	CHIRALCEL® OX-H	250	4.6	5	Analytical
63335	CHIRALCEL® OX-H	250	10	5	Semi-Prep
63337	CHIRALCEL® OX-H	20	10	5	Guard Semi-Prep
63345	CHIRALCEL® OX-H	250	20	5	Semi-Prep
63355	CHIRALCEL® OX-H	250	50	5	Preparative
63375	CHIRALCEL® OX-H	250	30	5	Semi-Prep
63394	CHIRALCEL® OX-H	150	2.1	5	Analytical
63423	CHIRALCEL® OX-H SFC	100	4.6	5	SFC Analytical
63435	CHIRALCEL® OX-H SFC	250	10	5	SFC Semi-Prep
63445	CHIRALCEL® OX-H SFC	250	20	5	SFC Semi-Prep
63455	CHIRALCEL® OX-H SFC	250	50	5	Preparative
63475	CHIRALCEL® OX-H SFC	250	30	5	SFC Semi-Prep

CHIRALCEL® OX-3

63511	CHIRALCEL® OX-3	10	4	3	Guard Cartridges (x 3)
63522	CHIRALCEL® OX-3	50	4.6	3	Analytical
63523	CHIRALCEL® OX-3	100	4.6	3	Analytical
63524	CHIRALCEL® OX-3	150	4.6	3	Analytical
63525	CHIRALCEL® OX-3	250	4.6	3	Analytical
63592	CHIRALCEL® OX-3	50	2.1	3	Analytical
63593	CHIRALCEL® OX-3	100	2.1	3	Analytical
63594	CHIRALCEL® OX-3	150	2.1	3	Analytical
63595	CHIRALCEL® OX-3	250	2.1	3	Analytical
63S82	CHIRALCEL® OX-3 SFC	50	3	3	SFC Analytical
63S83	CHIRALCEL® OX-3 SFC	100	3	3	SFC Analytical
63S84	CHIRALCEL® OX-3 SFC	150	3	3	SFC Analytical

CHIRALCEL® OX-RH

63711	CHIRALCEL® OX-RH	10	4	5	Guard Cartridges (x 3)
63724	CHIRALCEL® OX-RH	150	4.6	5	Analytical
63794	CHIRALCEL® OX-RH	150	2.1	5	Analytical

CHIRALCEL® OX-3R

63811	CHIRALCEL® OX-3R	10	4	3	Guard Cartridges (x 3)
63822	CHIRALCEL® OX-3R	50	4.6	3	Analytical
63823	CHIRALCEL® OX-3R	100	4.6	3	Analytical
63824	CHIRALCEL® OX-3R	150	4.6	3	Analytical

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
63892	CHIRALCEL® OX-3R	50	2.1	3	Analytical
63893	CHIRALCEL® OX-3R	100	2.1	3	Analytical
63894	CHIRALCEL® OX-3R	150	2.1	3	Analytical
63895	CHIRALCEL® OX-3R	250	2.1	3	Analytical

CHIRALCEL® OZ

42020	CHIRALCEL® OZ	-	Kg	20	CSP
42021	CHIRALCEL® OZ	-	100 g	20	CSP
42224	CHIRALCEL® OZ	150	4.6	20	Analytical
42225	CHIRALCEL® OZ	250	4.6	20	Analytical
4200E	CHIRALCEL® OZ-10	-	Kg	10	CSP
42G22	CHIRALCEL® OZ-10 SFC	50	4.6	10	SFC Analytical
42G24	CHIRALCEL® OZ-10 SFC	150	4.6	10	SFC Analytical
42G25	CHIRALCEL® OZ-10 SFC	250	4.6	10	SFC Analytical
42G94	CHIRALCEL® OZ-10 SFC	150	2.1	10	SFC Analytical
42G35	CHIRALCEL® OZ-10 SFC	250	10	10	SFC Semi-Prep
42G45	CHIRALCEL® OZ-10 SFC	250	20	10	SFC Semi-Prep
42G55	CHIRALCEL® OZ-10 SFC	250	50	10	SFC Preparative
42G75	CHIRALCEL® OZ-10 SFC	250	30	10	SFC Semi-Prep

CHIRALCEL® OZ-H

42311	CHIRALCEL® OZ-H	10	4	5	Guard Cartridges (x 3)
42321	CHIRALCEL® OZ-H	50	4.6	5	Analytical
42323	CHIRALCEL® OZ-H	100	4.6	5	Analytical
42324	CHIRALCEL® OZ-H	150	4.6	5	Analytical
42325	CHIRALCEL® OZ-H	250	4.6	5	Analytical
42335	CHIRALCEL® OZ-H	250	10	5	Semi-Prep
42337	CHIRALCEL® OZ-H	20	10	5	Guard Semi-Prep
42345	CHIRALCEL® OZ-H	250	20	5	Semi-Prep
42355	CHIRALCEL® OZ-H	250	50	5	Preparative
42375	CHIRALCEL® OZ-H	250	30	5	Semi-Prep
42394	CHIRALCEL® OZ-H	150	2.1	5	Analytical
42423	CHIRALCEL® OZ-H SFC	100	4.6	5	SFC Analytical
42435	CHIRALCEL® OZ-H SFC	250	10	5	SFC Semi-Prep
42445	CHIRALCEL® OZ-H SFC	250	20	5	SFC Semi-Prep
42455	CHIRALCEL® OZ-H SFC	250	50	5	SFC Preparative
42475	CHIRALCEL® OZ-H SFC	250	30	5	SFC Semi-Prep

CHIRALCEL® OZ-3

42511	CHIRALCEL® OZ-3	10	4	3	Guard Cartridges (x 3)
42522	CHIRALCEL® OZ-3	50	4.6	3	Analytical
42523	CHIRALCEL® OZ-3	100	4.6	3	Analytical
42524	CHIRALCEL® OZ-3	150	4.6	3	Analytical
42525	CHIRALCEL® OZ-3	250	4.6	3	Analytical
42592	CHIRALCEL® OZ-3	50	2.1	3	Analytical
42593	CHIRALCEL® OZ-3	100	2.1	3	Analytical
42594	CHIRALCEL® OZ-3	150	2.1	3	Analytical
42595	CHIRALCEL® OZ-3	250	2.1	3	Analytical
42S82	CHIRALCEL® OZ-3 SFC	50	3	3	SFC Analytical
42S83	CHIRALCEL® OZ-3 SFC	100	3	3	SFC Analytical
42S84	CHIRALCEL® OZ-3 SFC	150	3	3	SFC Analytical

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
42711	CHIRALCEL® OZ-RH	10	4	5	Guard Cartridges (x 3)
42724	CHIRALCEL® OZ-RH	150	4.6	5	Analytical
42794	CHIRALCEL® OZ-RH	150	2.1	5	Analytical

CHIRALCEL® OZ-3R

42811	CHIRALCEL® OZ-3R	10	4	3	Guard Cartridges (x 3)
42822	CHIRALCEL® OZ-3R	50	4.6	3	Analytical
42823	CHIRALCEL® OZ-3R	100	4.6	3	Analytical
42824	CHIRALCEL® OZ-3R	150	4.6	3	Analytical
42892	CHIRALCEL® OZ-3R	50	2.1	3	Analytical
42893	CHIRALCEL® OZ-3R	100	2.1	3	Analytical
42894	CHIRALCEL® OZ-3R	150	2.1	3	Analytical
42895	CHIRALCEL® OZ-3R	250	2.1	3	Analytical

CHIRALPAK® AD

19020	CHIRALPAK® AD	-	Kg	20	CSP
19021	CHIRALPAK® AD	-	100 g	20	CSP
19224	CHIRALPAK® AD	150	4.6	20	Analytical
19225	CHIRALPAK® AD	250	4.6	20	Analytical
19256	CHIRALPAK® AD	500	50	20	Preparative
19266	CHIRALPAK® AD	500	100	20	Preparative
19023	CHIRALPAK® AD	100	4.6	10	Analytical
19024	CHIRALPAK® AD	150	4.6	10	Analytical
19025	CHIRALPAK® AD	250	4.6	10	HPLC Analytical
19032	CHIRALPAK® AD	50	10	10	Guard Semi-Prep
19035	CHIRALPAK® AD	250	10	10	Semi-Prep
19042	CHIRALPAK® AD	50	20	10	Guard Semi-Prep
19045	CHIRALPAK® AD	250	20	10	Semi-Prep
19094	CHIRALPAK® AD	150	2.1	10	HPLC Analytical
1900E	CHIRALPAK® AD-10	-	Kg	10	CSP
19G22	CHIRALPAK® AD-10 SFC	50	4.6	10	SFC Analytical
19G24	CHIRALPAK® AD-10 SFC	150	4.6	10	SFC Analytical
19G25	CHIRALPAK® AD-10 SFC	250	4.6	10	SFC Analytical
19G94	CHIRALPAK® AD-10 SFC	150	2.1	10	SFC Analytical
19G35	CHIRALPAK® AD-10 SFC	250	10	10	SFC Semi-Prep
19G45	CHIRALPAK® AD-10 SFC	250	20	10	SFC Semi-Prep
19G55	CHIRALPAK® AD-10 SFC	250	50	10	SFC Preparative
19G75	CHIRALPAK® AD-10 SFC	250	30	10	SFC Semi-Prep

CHIRALPAK® AD-H

19311	CHIRALPAK® AD-H	10	4	5	Guard Cartridges (x 3)
19321	CHIRALPAK® AD-H	50	4.6	5	Analytical
19323	CHIRALPAK® AD-H	100	4.6	5	Analytical
19324	CHIRALPAK® AD-H	150	4.6	5	Analytical
19325	CHIRALPAK® AD-H	250	4.6	5	Analytical
19335	CHIRALPAK® AD-H	250	10	5	Semi-Prep
19337	CHIRALPAK® AD-H	20	10	5	Guard Semi-Prep
19345	CHIRALPAK® AD-H	250	20	5	Semi-Prep
19355	CHIRALPAK® AD-H	250	50	5	Preparative
19375	CHIRALPAK® AD-H	250	30	5	Semi-Prep
19394	CHIRALPAK® AD-H	150	2.1	5	Analytical

2018 Product List

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
19423	CHIRALPAK® AD-H SFC	100	4.6	5	SFC Analytical
19435	CHIRALPAK® AD-H SFC	250	10	5	SFC Semi-Prep
19445	CHIRALPAK® AD-H SFC	250	20	5	SFC Semi-Prep
19455	CHIRALPAK® AD-H SFC	250	50	5	SFC Preparative
19475	CHIRALPAK® AD-H SFC	250	30	5	SFC Semi-Prep

CHIRALPAK® AD-3

19511	CHIRALPAK® AD-3	10	4	3	Guard Cartridges (x 3)
19522	CHIRALPAK® AD-3	50	4.6	3	Analytical
19523	CHIRALPAK® AD-3	100	4.6	3	Analytical
19524	CHIRALPAK® AD-3	150	4.6	3	Analytical
19525	CHIRALPAK® AD-3	250	4.6	3	Analytical
19592	CHIRALPAK® AD-3	50	2.1	3	Analytical
19593	CHIRALPAK® AD-3	100	2.1	3	Analytical
19594	CHIRALPAK® AD-3	150	2.1	3	Analytical
19595	CHIRALPAK® AD-3	250	2.1	3	Analytical
19S82	CHIRALPAK® AD-3 SFC	50	3	3	SFC Analytical
19S83	CHIRALPAK® AD-3 SFC	100	3	3	SFC Analytical
19S84	CHIRALPAK® AD-3 SFC	150	3	3	SFC Analytical

CHIRALPAK® AD-RH

19711	CHIRALPAK® AD-RH	10	4	5	Guard Cartridges (x 3)
19724	CHIRALPAK® AD-RH	150	4.6	5	Analytical
19794	CHIRALPAK® AD-RH	150	2.1	5	Analytical

CHIRALPAK® AD-3R

19811	CHIRALPAK® AD-3R	10	4	3	Guard Cartridges (x 3)
19822	CHIRALPAK® AD-3R	50	4.6	3	Analytical
19823	CHIRALPAK® AD-3R	100	4.6	3	Analytical
19824	CHIRALPAK® AD-3R	150	4.6	3	Analytical
19892	CHIRALPAK® AD-3R	50	2.1	3	Analytical
19893	CHIRALPAK® AD-3R	100	2.1	3	Analytical
19894	CHIRALPAK® AD-3R	150	2.1	3	Analytical
19895	CHIRALPAK® AD-3R	250	2.1	3	Analytical

CHIRALPAK® AS-V

20020	CHIRALPAK® AS-V	-	Kg	20	CSP
20021	CHIRALPAK® AS-V	-	100 g	20	CSP
20224	CHIRALPAK® AS-V	150	4.6	20	Analytical
20225	CHIRALPAK® AS-V	250	4.6	20	Analytical
20256	CHIRALPAK® AS-V	500	50	20	Preparative
20266	CHIRALPAK® AS-V	500	100	20	Preparative

CHIRALPAK® AS

20024	CHIRALPAK® AS	150	4.6	10	Analytical
20025	CHIRALPAK® AS	250	4.6	10	HPLC Analytical
20032	CHIRALPAK® AS	50	10	10	Guard Semi-Prep
20035	CHIRALPAK® AS	250	10	10	Semi-Prep
20042	CHIRALPAK® AS	50	20	10	Guard Semi-Prep
20045	CHIRALPAK® AS	250	20	10	Semi-Prep
20094	CHIRALPAK® AS	150	2.1	10	Analytical
2000E	CHIRALPAK® AS-10	-	Kg	10	CSP

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
20G22	CHIRALPAK® AS-10 SFC	50	4.6	10	SFC Analytical
20G24	CHIRALPAK® AS-10 SFC	150	4.6	10	SFC Analytical
20G25	CHIRALPAK® AS-10 SFC	250	4.6	10	SFC Analytical
20G94	CHIRALPAK® AS-10 SFC	150	2.1	10	SFC Analytical
20G35	CHIRALPAK® AS-10 SFC	250	10	10	SFC Semi-Prep
20G45	CHIRALPAK® AS-10 SFC	250	20	10	SFC Semi-Prep
20G55	CHIRALPAK® AS-10 SFC	250	50	10	SFC Preparative
20G75	CHIRALPAK® AS-10 SFC	250	30	10	SFC Semi-Prep

CHIRALPAK® AS-H

20311	CHIRALPAK® AS-H	10	4	5	Guard Cartridges (x 3)
20321	CHIRALPAK® AS-H	50	4.6	5	Analytical
20323	CHIRALPAK® AS-H	100	4.6	5	Analytical
20324	CHIRALPAK® AS-H	150	4.6	5	Analytical
20325	CHIRALPAK® AS-H	250	4.6	5	Analytical
20335	CHIRALPAK® AS-H	250	10	5	Semi-Prep
20337	CHIRALPAK® AS-H	20	10	5	Guard Semi-Prep
20345	CHIRALPAK® AS-H	250	20	5	Semi-Prep
20355	CHIRALPAK® AS-H	250	50	5	Preparative
20375	CHIRALPAK® AS-H	250	30	5	Semi-Prep
20394	CHIRALPAK® AS-H	150	2.1	5	Analytical
20423	CHIRALPAK® AS-H SFC	100	4.6	5	SFC Analytical
20435	CHIRALPAK® AS-H SFC	250	10	5	SFC Semi-Prep
20445	CHIRALPAK® AS-H SFC	250	20	5	SFC Semi-Prep
20455	CHIRALPAK® AS-H SFC	250	50	5	SFC Preparative
20475	CHIRALPAK® AS-H SFC	250	30	5	SFC Semi-Prep

CHIRALPAK® AS-3

20511	CHIRALPAK® AS-3	10	4	3	Guard Cartridges (x 3)
20522	CHIRALPAK® AS-3	50	4.6	3	Analytical
20523	CHIRALPAK® AS-3	100	4.6	3	Analytical
20524	CHIRALPAK® AS-3	150	4.6	3	Analytical
20525	CHIRALPAK® AS-3	250	4.6	3	Analytical
20592	CHIRALPAK® AS-3	50	2.1	3	Analytical
20593	CHIRALPAK® AS-3	100	2.1	3	Analytical
20594	CHIRALPAK® AS-3	150	2.1	3	Analytical
20595	CHIRALPAK® AS-3	250	2.1	3	Analytical
20S82	CHIRALPAK® AS-3 SFC	50	3	3	SFC Analytical
20S83	CHIRALPAK® AS-3 SFC	100	3	3	SFC Analytical
20S84	CHIRALPAK® AS-3 SFC	150	3	3	SFC Analytical

CHIRALPAK® AS-RH

20711	CHIRALPAK® AS-RH	10	4	5	Guard Cartridges (x 3)
20724	CHIRALPAK® AS-RH	150	4.6	5	Analytical
20794	CHIRALPAK® AS-RH	150	2.1	5	Analytical

CHIRALPAK® AS-3R

20811	CHIRALPAK® AS-3R	10	4	3	Guard Cartridges (x 3)
20822	CHIRALPAK® AS-3R	50	4.6	3	Analytical
20823	CHIRALPAK® AS-3R	100	4.6	3	Analytical
20824	CHIRALPAK® AS-3R	150	4.6	3	Analytical

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
20892	CHIRALPAK® AS-3R	50	2.1	3	Analytical
20893	CHIRALPAK® AS-3R	100	2.1	3	Analytical
20894	CHIRALPAK® AS-3R	150	2.1	3	Analytical
20895	CHIRALPAK® AS-3R	250	2.1	3	Analytical

CHIRALPAK® AY

47020	CHIRALPAK® AY	-	Kg	20	CSP
47021	CHIRALPAK® AY	-	100 g	20	CSP
47224	CHIRALPAK® AY	150	4.6	20	Analytical
47225	CHIRALPAK® AY	250	4.6	20	Analytical
47024	CHIRALPAK® AY	150	4.6	10	Analytical
47025	CHIRALPAK® AY	250	4.6	10	Analytical
47135	CHIRALPAK® AY SFC	250	10	10	SFC Semi-Prep
47145	CHIRALPAK® AY SFC	250	20	10	SFC Semi-Prep
47175	CHIRALPAK® AY SFC	250	30	10	SFC Semi-Prep
4700E	CHIRALPAK® AY-10	-	Kg	10	CSP
47G22	CHIRALPAK® AY-10 SFC	50	4.6	10	SFC Analytical
47G24	CHIRALPAK® AY-10 SFC	150	4.6	10	SFC Analytical
47G25	CHIRALPAK® AY-10 SFC	250	4.6	10	SFC Analytical
47G94	CHIRALPAK® AY-10 SFC	150	2.1	10	SFC Analytical
47G35	CHIRALPAK® AY-10 SFC	250	10	10	SFC Semi-Prep
47G45	CHIRALPAK® AY-10 SFC	250	20	10	SFC Semi-Prep
47G55	CHIRALPAK® AY-10 SFC	250	50	10	SFC Preparative
47G75	CHIRALPAK® AY-10 SFC	250	30	10	SFC Semi-Prep

CHIRALPAK® AY-H

47311	CHIRALPAK® AY-H	10	4	5	Guard Cartridges (x 3)
47321	CHIRALPAK® AY-H	50	4.6	5	Analytical
47323	CHIRALPAK® AY-H	100	4.6	5	Analytical
47324	CHIRALPAK® AY-H	150	4.6	5	Analytical
47325	CHIRALPAK® AY-H	250	4.6	5	Analytical
47335	CHIRALPAK® AY-H	250	10	5	Semi-Prep
47337	CHIRALPAK® AY-H	20	10	5	Guard Semi-Prep
47345	CHIRALPAK® AY-H	250	20	5	Semi-Prep
47355	CHIRALPAK® AY-H	250	50	5	Preparative
47375	CHIRALPAK® AY-H	250	30	5	Semi-Prep
47394	CHIRALPAK® AY-H	150	2.1	5	Analytical
47423	CHIRALPAK® AY-H SFC	100	4.6	5	SFC Analytical
47435	CHIRALPAK® AY-H SFC	250	10	5	SFC Semi-Prep
47445	CHIRALPAK® AY-H SFC	250	20	5	SFC Semi-Prep
47475	CHIRALPAK® AY-H SFC	250	30	5	SFC Semi-Prep
47455	CHIRALPAK® AY-H SFC	250	50	5	SFC Preparative

CHIRALPAK® AY-3

47511	CHIRALPAK® AY-3	10	4	3	Guard Cartridges (x 3)
47522	CHIRALPAK® AY-3	50	4.6	3	Analytical
47523	CHIRALPAK® AY-3	100	4.6	3	Analytical
47524	CHIRALPAK® AY-3	150	4.6	3	Analytical
47525	CHIRALPAK® AY-3	250	4.6	3	Analytical
47592	CHIRALPAK® AY-3	50	2.1	3	Analytical
47593	CHIRALPAK® AY-3	100	2.1	3	Analytical
47594	CHIRALPAK® AY-3	150	2.1	3	Analytical
47594	CHIRALPAK® AY-3	250	2.1	3	Analytical

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
47S82	CHIRALPAK® AY-3 SFC	50	3	3	SFC Analytical
47S83	CHIRALPAK® AY-3 SFC	100	3	3	SFC Analytical
47S84	CHIRALPAK® AY-3 SFC	150	3	3	SFC Analytical

CHIRALPAK® AY-RH

47711	CHIRALPAK® AY-RH	10	4	5	Guard Cartridges (x 3)
47724	CHIRALPAK® AY-RH	150	4.6	5	Analytical
47794	CHIRALPAK® AY-RH	150	2.1	5	Analytical

CHIRALPAK® AY-3R

47811	CHIRALPAK® AY-3R	10	4	3	Guard Cartridges (x 3)
47822	CHIRALPAK® AY-3R	50	4.6	3	Analytical
47823	CHIRALPAK® AY-3R	100	4.6	3	Analytical
47824	CHIRALPAK® AY-3R	150	4.6	3	Analytical
47892	CHIRALPAK® AY-3R	50	2.1	3	Analytical
47893	CHIRALPAK® AY-3R	100	2.1	3	Analytical
47894	CHIRALPAK® AY-3R	150	2.1	3	Analytical
47895	CHIRALPAK® AY-3R	250	2.1	3	Analytical

CHIRALPAK® AZ

61020	CHIRALPAK® AZ	-	Kg	20	CSP
61021	CHIRALPAK® AZ	-	100 g	20	CSP
61224	CHIRALPAK® AZ	150	4.6	20	Analytical
61225	CHIRALPAK® AZ	250	4.6	20	Analytical
61024	CHIRALPAK® AZ	150	4.6	10	Analytical
61025	CHIRALPAK® AZ	250	4.6	10	Analytical
61135	CHIRALPAK® AZ	250	10	10	Semi-Prep
61145	CHIRALPAK® AZ	250	20	10	Semi-Prep
61175	CHIRALPAK® AZ	250	30	10	Semi-Prep
6100E	CHIRALPAK® AZ-10	-	Kg	10	CSP
61G22	CHIRALPAK® AZ-10 SFC	50	4.6	10	SFC Analytical
61G24	CHIRALPAK® AZ-10 SFC	150	4.6	10	SFC Analytical
61G25	CHIRALPAK® AZ-10 SFC	250	4.6	10	SFC Analytical
61G94	CHIRALPAK® AZ-10 SFC	150	2.1	10	SFC Analytical
61G35	CHIRALPAK® AZ-10 SFC	250	10	10	SFC Semi-Prep
61G45	CHIRALPAK® AZ-10 SFC	250	20	10	SFC Semi-Prep
61G55	CHIRALPAK® AZ-10 SFC	250	50	10	SFC Preparative
61G75	CHIRALPAK® AZ-10 SFC	250	30	10	SFC Semi-Prep

CHIRALPAK® AZ-H

61311	CHIRALPAK® AZ-H	10	4	5	Guard Cartridges (x 3)
61321	CHIRALPAK® AZ-H	50	4.6	5	Analytical
61323	CHIRALPAK® AZ-H	100	4.6	5	Analytical
61324	CHIRALPAK® AZ-H	150	4.6	5	Analytical
61325	CHIRALPAK® AZ-H	250	4.6	5	Analytical
61335	CHIRALPAK® AZ-H	250	10	5	Semi-Prep
61337	CHIRALPAK® AZ-H	20	10	5	Guard Semi-Prep
61345	CHIRALPAK® AZ-H	250	20	5	Semi-Prep
61355	CHIRALPAK® AZ-H	250	50	5	Preparative
61375	CHIRALPAK® AZ-H	250	30	5	Semi-Prep
61394	CHIRALPAK® AZ-H	150	2.1	5	Analytical
61423	CHIRALPAK® AZ-H SFC	100	4.6	5	SFC Analytical
61435	CHIRALPAK® AZ-H SFC	250	10	5	SFC Semi-Prep

2018 Product List

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
61445	CHIRALPAK® AZ-H SFC	250	20	5	SFC Semi-Prep
61455	CHIRALPAK® AZ-H SFC	250	50	5	SFC Preparative
61475	CHIRALPAK® AZ-H SFC	250	30	5	SFC Semi-Prep

CHIRALPAK® AZ-3

61511	CHIRALPAK® AZ-3	10	4	3	Guard Cartridges (x 3)
61522	CHIRALPAK® AZ-3	50	4.6	3	Analytical
61523	CHIRALPAK® AZ-3	100	4.6	3	Analytical
61524	CHIRALPAK® AZ-3	150	4.6	3	Analytical
61525	CHIRALPAK® AZ-3	250	4.6	3	Analytical
61592	CHIRALPAK® AZ-3	50	2.1	3	Analytical
61593	CHIRALPAK® AZ-3	100	2.1	3	Analytical
61594	CHIRALPAK® AZ-3	150	2.1	3	Analytical
61595	CHIRALPAK® AZ-3	250	2.1	3	Analytical
61S82	CHIRALPAK® AZ-3 SFC	50	3	3	SFC Analytical
61S83	CHIRALPAK® AZ-3 SFC	100	3	3	SFC Analytical
61S84	CHIRALPAK® AZ-3 SFC	150	3	3	SFC Analytical

CHIRALPAK® AZ-RH

61711	CHIRALPAK® AZ-RH	10	4	5	Guard Cartridges (x 3)
61724	CHIRALPAK® AZ-RH	150	4.6	5	Analytical
61794	CHIRALPAK® AZ-RH	150	2.1	5	Analytical

CHIRALPAK® AZ-3R

61811	CHIRALPAK® AZ-3R	10	4	3	Guard Cartridges (x 3)
61822	CHIRALPAK® AZ-3R	50	4.6	3	Analytical
61823	CHIRALPAK® AZ-3R	100	4.6	3	Analytical
61824	CHIRALPAK® AZ-3R	150	4.6	3	Analytical
61892	CHIRALPAK® AZ-3R	50	2.1	3	Analytical
61893	CHIRALPAK® AZ-3R	100	2.1	3	Analytical
61894	CHIRALPAK® AZ-3R	150	2.1	3	Analytical
61895	CHIRALPAK® AZ-3R	250	2.1	3	Analytical

MPLC COLUMNS - METHOD DEVELOPMENT

80223	CHIRALPAK® IA	100	4.6	20	Method Development
83223	CHIRALPAK® IC	100	4.6	20	Method Development
84223	CHIRALPAK® ID	100	4.6	20	Method Development
85223	CHIRALPAK® IE	100	4.6	20	Method Development
86223	CHIRALPAK® IF	100	4.6	20	Method Development

MPLC COLUMNS

80M73	CHIRALFLASH IA	100	30	20	Preparative
83M73	CHIRALFLASH IC	100	30	20	Preparative
84M73	CHIRALFLASH ID	100	30	20	Preparative
85M73	CHIRALFLASH IE	100	30	20	Preparative
86M73	CHIRALFLASH IF	100	30	20	Preparative

CHIRAL TLC PLATES

5ST5A	2L-ChiralTLC (5-pkg/1 each)	100	200	20	Method Development
80T2A	2L-ChiralTLC IA (2-pkg)	100	200	20	Method Development
83T2A	2L-ChiralTLC IC (2-pkg)	100	200	20	Method Development
84T2A	2L-ChiralTLC ID (2-pkg)	100	200	20	Method Development

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
85T2A	2L-ChiralTLC IE (2-pkg)	100	200	20	Method Development
86T2A	2L-ChiralTLC IF (2-pkg)	100	200	20	Method Development

PROTEIN-BASED COLUMNS

CHIRALPAK® AGP

30713	CHIRALPAK® AGP	100	4	5	Analytical
30714	CHIRALPAK® AGP	150	4	5	Analytical
30712	CHIRALPAK® AGP	50	4	5	Analytical
30783	CHIRALPAK® AGP	100	3	5	Analytical
30784	CHIRALPAK® AGP	150	3	5	Analytical
30782	CHIRALPAK® AGP	50	3	5	Analytical
30793	CHIRALPAK® AGP	100	2	5	Analytical
30794	CHIRALPAK® AGP	150	2	5	Analytical
30792	CHIRALPAK® AGP	50	2	5	Analytical
30733	CHIRALPAK® AGP	100	10	5	Semi-Prep
30734	CHIRALPAK® AGP	150	10	5	Semi-Prep
30781	CHIRALPAK® AGP (2-pkg)	10	3	5	Guard Cartridges (x 2)
30791	CHIRALPAK® AGP (2-pkg)	10	2	5	Guard Cartridges (x 2)
30711	CHIRALPAK® AGP (2-pkg)	10	4	5	Guard Cartridges (x 2)

CHIRALPAK® CBH

33713	CHIRALPAK® CBH	100	4	5	Analytical
33714	CHIRALPAK® CBH	150	4	5	Analytical
33712	CHIRALPAK® CBH	50	4	5	Analytical
33783	CHIRALPAK® CBH	100	3	5	Analytical
33784	CHIRALPAK® CBH	150	3	5	Analytical
33782	CHIRALPAK® CBH	50	3	5	Analytical
33793	CHIRALPAK® CBH	100	2	5	Analytical
33794	CHIRALPAK® CBH	150	2	5	Analytical
33792	CHIRALPAK® CBH	50	2	5	Analytical
33733	CHIRALPAK® CBH	100	10	5	Semi-Prep
33734	CHIRALPAK® CBH	150	10	5	Semi-Prep
33781	CHIRALPAK® CBH (2-pkg)	10	3	5	Guard Cartridges (x 2)
33791	CHIRALPAK® CBH (2-pkg)	10	2	5	Guard Cartridges (x 2)
33711	CHIRALPAK® CBH (2-pkg)	10	4	5	Guard Cartridges (x 2)

CHIRALPAK® HSA

34713	CHIRALPAK® HSA	100	4	5	Analytical
34714	CHIRALPAK® HSA	150	4	5	Analytical
34712	CHIRALPAK® HSA	50	4	5	Analytical
34783	CHIRALPAK® HSA	100	3	5	Analytical
34784	CHIRALPAK® HSA	150	3	5	Analytical
34782	CHIRALPAK® HSA	50	3	5	Analytical
34793	CHIRALPAK® HSA	100	2	5	Analytical
34794	CHIRALPAK® HSA	150	2	5	Analytical
34792	CHIRALPAK® HSA	50	2	5	Analytical
34733	CHIRALPAK® HSA	100	10	5	Semi-Prep
34734	CHIRALPAK® HSA	150	10	5	Semi-Prep
34781	CHIRALPAK® HSA (2-pkg)	10	3	5	Guard Cartridges (x 2)
34791	CHIRALPAK® HSA (2-pkg)	10	2	5	Guard Cartridges (x 2)
34711	CHIRALPAK® HSA (2-pkg)	10	4	5	Guard Cartridges (x 2)

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
CHIRALPAK® BSA					
35712	CHIRALPAK® BSA	50	4	5	Analytical
35714	CHIRALPAK® BSA	150	4	5	Analytical
35782	CHIRALPAK® BSA	50	3	5	Analytical
35792	CHIRALPAK® BSA	50	2	5	Analytical
35781	CHIRALPAK® BSA	10	3	5	Guard Cartridges (x 2)
35791	CHIRALPAK® BSA	10	2	5	Guard Cartridges (x 2)
35711	CHIRALPAK® BSA	10	4	5	Guard Cartridges (x 2)

CHIRALPAK® DSA

36712	CHIRALPAK® DSA	50	4	5	Analytical
36714	CHIRALPAK® DSA	150	4	5	Analytical
36782	CHIRALPAK® DSA	50	3	5	Analytical
36792	CHIRALPAK® DSA	50	2	5	Analytical
36781	CHIRALPAK® DSA	10	3	5	Guard Cartridges (x 2)
36791	CHIRALPAK® DSA	10	2	5	Guard Cartridges (x 2)
36711	CHIRALPAK® DSA	10	4	5	Guard Cartridges (x 2)

CHIRALPAK® MSA

37712	CHIRALPAK® MSA	50	4	5	Analytical
37782	CHIRALPAK® MSA	50	3	5	Analytical

CHIRALPAK® PSA

43712	CHIRALPAK® PSA	50	4	5	Analytical
43782	CHIRALPAK® PSA	50	3	5	Analytical

CHIRALPAK® RSA

38712	CHIRALPAK® RSA	50	4	5	Analytical
38714	CHIRALPAK® RSA	150	4	5	Analytical
38782	CHIRALPAK® RSA	50	3	5	Analytical
38792	CHIRALPAK® RSA	50	2	5	Analytical
38781	CHIRALPAK® RSA	10	3	5	Guard Cartridges (x 2)
38791	CHIRALPAK® RSA	10	2	5	Guard Cartridges (x 2)
38711	CHIRALPAK® RSA	10	4	5	Guard Cartridges (x 2)

PROTEIN-BASED COLUMN ACCESSORIES

00081	GUARD CARTRIDGE HOLDER*	Holder for Guard Cartridges			
000D1	GUARD COLUMN COUPLER	Coupler			
000D2	GUARD COLUMN COUPLER MICRO	Coupler			

* for protein-based guard cartridges

CHIRAL SELECTORS FOR SPECIAL APPLICATIONS**IMMOBILISED CROWNPAK® CR-I(+) & (-)**

53784	CROWNPAK® CR-I(+)	150	3	5	Analytical
54784	CROWNPAK® CR-I(-)	150	3	5	Analytical

COATED CROWNPAK® CR(+) & (-)

27711	CROWNPAK® CR	10	4	5	Guard
27714	CROWNPAK® CR(+)	150	4	5	Analytical
28714	CROWNPAK® CR(-)	150	4	5	Analytical

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
CHIRALPAK® ZWIX(+) & ZWIX(-)					
51334	CHIRALPAK® ZWIX(+)	150	10	5	Semi-Prep
51335	CHIRALPAK® ZWIX(+)	250	10	5	Semi-Prep
51344	CHIRALPAK® ZWIX(+)	150	21	5	Semi-Prep
51345	CHIRALPAK® ZWIX(+)	250	21	5	Semi-Prep
51514	CHIRALPAK® ZWIX(+)	150	4	3	Analytical
51515	CHIRALPAK® ZWIX(+)	250	4	3	Analytical
51584	CHIRALPAK® ZWIX(+)	150	3	3	Analytical
51585	CHIRALPAK® ZWIX(+)	250	3	3	Analytical
52334	CHIRALPAK® ZWIX(-)	150	10	5	Semi-Prep
52335	CHIRALPAK® ZWIX(-)	250	10	5	Semi-Prep
52344	CHIRALPAK® ZWIX(-)	150	21	5	Semi-Prep
52345	CHIRALPAK® ZWIX(-)	250	21	5	Semi-Prep
52514	CHIRALPAK® ZWIX(-)	150	4	3	Analytical
52515	CHIRALPAK® ZWIX(-)	250	4	3	Analytical
52584	CHIRALPAK® ZWIX(-)	150	3	3	Analytical
52585	CHIRALPAK® ZWIX(-)	250	3	3	Analytical

CHIRALPAK® QD-AX & QN-AX

31311	CHIRALPAK® QD-AX	10	4	5	Guard Cartridges (x 3)
31324	CHIRALPAK® QD-AX	150	4.6	5	Analytical
31394	CHIRALPAK® QD-AX	150	2.1	5	Analytical
31344	CHIRALPAK® QD-AX	150	20	5	Semi-Prep
31444	CHIRALPAK® QD-AX SFC	150	20	5	SFC Semi-Prep
32311	CHIRALPAK® QN-AX	10	4	5	Guard Cartridges (x 3)
32324	CHIRALPAK® QN-AX	150	4.6	5	Analytical
32394	CHIRALPAK® QN-AX	150	2.1	5	Analytical
32344	CHIRALPAK® QN-AX	150	20	5	Semi-Prep
32444	CHIRALPAK® QN-AX SFC	150	20	5	SFC Semi-Prep

CHIRALPAK® MA(+)

21822	CHIRALPAK® MA(+)	50	4.6	3	Analytical
-------	------------------	----	-----	---	------------

CHIRALPAK® WH

25622	CHIRALPAK® WH	50	4.6	10	Guard
25625	CHIRALPAK® WH	250	4.6	10	Analytical

SFC ACHIRAL COLUMNS**DAICEL DCPAK® SFC-A**

AA424	DAICEL DCPak® SFC-A	150	4.6	5	SFC Analytical
AA425	DAICEL DCPak® SFC-A	250	4.6	5	SFC Analytical
AA435	DAICEL DCPak® SFC-A	250	10	5	SFC Semi-Prep
AA445	DAICEL DCPak® SFC-A	250	21	5	SFC Semi-Prep
AA494	DAICEL DCPak® SFC-A	150	2.1	5	SFC Analytical
AA495	DAICEL DCPak® SFC-A	250	2.1	5	SFC Analytical
AAS82	DAICEL DCPak® SFC-A	50	3	3	SFC Analytical
AAS83	DAICEL DCPak® SFC-A	100	3	3	SFC Analytical
AAS84	DAICEL DCPak® SFC-A	150	3	3	SFC Analytical
AAS92	DAICEL DCPak® SFC-A	50	2.1	3	SFC Analytical
AAS93	DAICEL DCPak® SFC-A	100	2.1	3	SFC Analytical
AAS94	DAICEL DCPak® SFC-A	150	2.1	3	SFC Analytical

2018 Product List

REFERENCE	PHASE	LENGTH (mm)	ID (mm)	PARTICLE SIZE (μ)	TYPE
DAICEL DCpak® SFC-B					
AB424	DAICEL DCpak® SFC-B	150	4.6	5	SFC Analytical
AB425	DAICEL DCpak® SFC-B	250	4.6	5	SFC Analytical
AB435	DAICEL DCpak® SFC-B	250	10	5	SFC Semi-Prep
AB445	DAICEL DCpak® SFC-B	250	21	5	SFC Semi-Prep
AB494	DAICEL DCpak® SFC-B	150	2.1	5	SFC Analytical
AB495	DAICEL DCpak® SFC-B	250	2.1	5	SFC Analytical
ABS82	DAICEL DCpak® SFC-B	50	3	3	SFC Analytical
ABS83	DAICEL DCpak® SFC-B	100	3	3	SFC Analytical
ABS84	DAICEL DCpak® SFC-B	150	3	3	SFC Analytical
ABS92	DAICEL DCpak® SFC-B	50	2.1	3	SFC Analytical
ABS93	DAICEL DCpak® SFC-B	100	2.1	3	SFC Analytical
ABS94	DAICEL DCpak® SFC-B	150	2.1	3	SFC Analytical

COLUMN ACCESSORIES

00011	GUARD CARTRIDGE HOLDER†	-	-	-	Holder for Guard Cartridges
00024	Column Jacket	150	4.6	-	Hardware
00025	Column Jacket	250	4.6	-	Hardware
000D1- DAIC	Guard Column Coupler	-	-	-	Hardware
FK002	Guard Filter Kit (1 holder, 3 glass filters, 2 fittings and 1 tubing)	-	-	-	Guard Filter
FR002	Replacement Guard Filters (x 6)	-	-	-	Guard Filters

†For all guard cartridges



CHIRAL TECHNOLOGIES EUROPE SAS

Parc d'Innovation

160, Bd Gonthier d'Andernach - CS 80140
67404 Illkirch Cedex FRANCE

Tel: +33 (0) 3 88 79 52 00

Fax: +33 (0) 3 88 66 71 66

Email: cte@chiral.fr

WWW.CHIRALTECH.COM

CHIRALPAK, CHIRALCEL,
and CROWNPAK are registered
trademarks of DAICEL CORPORATION.

© 2018 CHIRAL TECHNOLOGIES EUROPE